
OPTICAL MODULE ASSEMBLY INSTRUCTIONS



Andrey Chirikov, Brady Kreil, Gus Hawley, Harry Zhou, Mikhail Ignatyev, Stephanie Monagan

UBC Manufacturing Engineering Capstone Project 2024/2025

Optical Module Assembly Instructions Rev. 1, 09-Apr-2025

1.	Introduction	1	
2.	Design Overview	1	
3.	BOM	2	
4.	Machined Parts	2	
	4.1 Frame	2	
	4.2 Plates	3	,
5.	3D Printed Parts	4	
	5.1 Overview	4	
	5.2 Recommended Settings	5)
	5.3 Camera Carriage	5	,
	5.4 Mirrored Brackets	6)
6.	Assembly	6	į
	6.1 Frame Assembly	6	,
	6.2 Camera Carriage	7	
	6.3 Belt Assembly	8	,
	6.4 Platform Assembly	10	l
	6.5 Linear Rail Alignment	11	
	6.6 Z Motor and Ball Screw Alignment	12	
	6.7 Electrical Assembly	12	
	6.8 Overview	12	
	6.9 E-STOP	13	,
	6.10 Buck Converters	13	,
7.	Software Setup	14	
	7.1 Raspberry Pi Settings	14	
	7.2 Networking	14	
	7.3 Windows PC	16	,
	7.3.1 Ethernet	16	,
	7.3.2 Folder Directories	16	,
	7.3.3 Fiji (ImageJ)	16	į
Αį	ppendix	19	,
	A.1 – Full Parts List	19	1
	A.2 – List of Fasteners per Part	23	į

Optical Module Assembly Instructions Rev. 1, 09-Apr-2025

A.3 – List of Machined Parts	26
A.4 – Assembly List - Electronics	27
A.5 – BOM - Electronics	29

1. Introduction

This document is intended to go over how to assemble the optical module from scratch. The optical module is an automated microscope, designed to be integrated in a self-driving lab. Many components and information will be stored on our GitHub and referenced throughout the document:

https://github.com/ManuCapstone1/OpticalModule/tree/main

2. Design Overview

This section provides an overview of the final design. The entire assembly is mounted on 20x20 aluminum extrusions, forming a rigid structural frame. The sample platform is a breadboard, allowing flexible placement of the sample at any position. The platform assembly is supported by two trolleys running along rails, enhancing the rigidity of the structure. It is connected to a ball screw mechanism actuated by a DC motor to ensure precise movement along the Z-axis.

The optical module moves in the XY plane and is mounted on a trolley rail system, similar to the Z-axis setup. However, its movement is controlled by a CoreXY motion system, which uses belts and DC motors for efficient and precise positioning. The CAD model shown in Figure 1 shows the complete design, excluding belts and fasteners. A detailed list of fasteners and their connection points is provided in Appendix A.2.






Figure 1: Front and back view of the main assembly of the project

3. BOM

The full list of all the needed parts and tools are found in appendix A.1 – Full Parts ListA.1 – Full . Parts lists for fasteners and machined parts are also in appendix A. For a more detailed and organized view, an excel spreadsheet is in the GitHub repository in OpticalModule/Documentation/Github BOM.

4. Machined Parts

4.1 Frame

The extrusions used for the frame of the module need to be extremely precise and have a flat surface on each end to ensure a square fit between all of the extrusions.

The recommended cutting procedure is as follows:

- 1. Measure the required extrusion length, adding 5 mm to the specified length. Mark the measurement on the extrusion using a marker. Measure from the factory cut end, when possible, to ensure better cut surface quality.
- 2. Cut the extrusion using a bandsaw.
- 3. Set up a milling machine for cutting the extrusions. Secure an adjustable end stop to serve as a reference point for all extrusion pieces. Keep the end stop fixed for pieces of the same length, but adjust it as needed for other lengths.

- 4. Place the extrusion in a vise, press it against the end stop, and tighten the vise securely.
- 5. Load an edge finder tool, turn on the spindle, and carefully move the tool against the extrusion's end-stop edge from the side until it shifts—this marks the zero point, factoring in the edge finder's radius.
- 6. Raise the edge finder above the extrusion piece's edge and verify its position. It should match the radius of the edge finder. If correct, zero the coordinates again—this is the final zero point.
- 7. Replace the edge finder with a milling tool. For 20x20 extrusions, a 1" tool is recommended. Calculate the final milling position as the desired length plus the tool's radius.
- 8. Mill the free end of the extrusion using 0.5 mm passes, finishing with a down milling cut for a better surface finish. Verify the final length, ensuring it falls within the 0.1 mm tolerance.

4.2 Plates

The recommended machining procedures for the enclosure backplate, platform baseplate and the Breadboard Platform are described below. For dimensions and other information, reference Appendix A.3 – List of Machined Parts.

Enclosure Backplate

- 1. Select a metal sheet of the specified thickness. Using a marker, outline the specified length and width according to the drawing.
- 2. Cut the sheet to size using a shear cutter. For better accuracy, cut slightly beyond the marked dimensions to allow for final adjustments by grinding if needed.
- 3. Place the sheet on a flat granite surface and use a height gauge to scribe lines indicating the specified hole positions.
- 4. At each intersection of the scribed lines, use a center punch and hammer to create small dents. This will help guide the drill bit accurately.
- 5. Prepare M2.5, M3, and M5 drill bits for use with a drill press. Secure the sheet face-up in a vise, ensuring it is tightly clamped.

- 6. Load the first drill bit into the drill press. Align it with a punched dent and lower it slightly to confirm proper positioning. Adjust the vise if necessary.
 - a. Drill through each marked hole, following the specified locations.
 - b. Repeat step 6 for each drill bit size.
- 7. Use a deburring tool to remove any burrs from the back side of the plate.

Base plate

The base plate was fabricated using a waterjet cutter. A correctly sized metal plate was positioned in the machine, and the DXF file containing the final shape was uploaded to the cutter's computer. The plate was then cut according to standard waterjet cutting procedures.

Breadboard Platform

- 1. Place the original Breadboard Platform on a flat granite surface. Use a height gauge to scribe lines marking the specified hole positions.
- 2. At each intersection of the scribed lines, use a center punch and hammer to create small dents to guide the drill bit.
- 3. Secure the Breadboard Plate firmly in a vise.
- 4. Select an M6 drill bit suitable for use with a drill press. Load it into the drill press and lower it onto a dent to verify proper alignment. Adjust the vise if necessary.
- 5. Drill through each marked hole.
- 6. Use a deburring tool to remove burrs from the opposite side of the plate.

5. 3D Printed Parts

5.1 Overview

All the .STL files for the parts that need to be printed can be found on the GitHub under OpticalModule/CAD/STL Files/. To see a list of the parts go to the "3D Printed" category in appendix A.1 – Full Parts List. Furthermore, a list of all the fasteners corresponding to each part is in appendix

A.2 – List of Fasteners per Part.

5.2 Recommended Settings

Below are the settings and parameters the team used in the slicer to print the 3D printed parts with a FDM printer. The parts were printed with the intention of having good structural integrity and rigidity to prolong the lifespan and reduce vibrations. These settings are *suggestions*, and users are at liberty to modify them as needed.

One exception to these settings are the jigs for the linear rails. Since they are not used in the final product, nor structurally, standard/weaker settings can be used instead.

Parameter	Setting/Parameter	Notes
Material	ABS, PLA, PLA+, PETG	Stronger materials will work as well but are not needed.
Resolution	0.02 mm	
Infill	~50%	
Infill Pattern	Cubic	Good for strength in all directions.
Supports	>45 degrees	
Shell Thickness	1.2 mm in both directions.	

Table 1: Recommended printer settings for 3D printed parts.

5.3 Camera Carriage

The STL file "EnclosureBase" is the main body of the camera carriage. Orient the part as shown below in Figure 1gure 2, as it ensures a nice outer finish and minimal shifting errors and cracking. When the print is complete, removing the supports is difficult. If needed, use pliers, a Dremel and/or similar tools.




Figure 1: Orientation of how the "EnclosureBase" STL file should be placed on the slicer.

5.4 Mirrored Brackets

Two STL files "XY Motor Mount" and "Z Carriage" need two parts each, one as is, and the other mirrored. Use the slicer to mirror one of the parts. Figure 32 and Figure 34 show how the parts should be mirrored. Note that the parts in the figures were placed to show how the parts should be mirrored, *not* how the parts should be oriented for printing.




Figure 32: View of how "XY Motor Mount" STL files should be mirrored in the slicer.




Figure 3: View of how "Z Carriage" STL File should be mirrored in the slicer.

6. Assembly

6.1 Frame Assembly

- To assemble the frame, start by constructing the sides. Lay two 465 mm extrusions on a flat surface and position three 370 mm extrusions perpendicularly between them.
- 2. Secure the 370 mm pieces to one of the 465 mm extrusions using M5x10 fasteners and M5 T-nuts (detailed list of fasteners in Appendix A.2), then repeat for the other 465 mm extrusion.
- Check squareness by standing the subassembly vertically and using a bubble gauge to verify perpendicularity. Adjust the brackets as needed to center the bubble.

- 4. Next, lay assembled side flat and attach four 370 mm extrusions vertically using M5x10 fasteners and M5 T-nuts, ensuring proper alignment.
- 5. Attach four brackets—two on each side—to the 395 mm extrusion, then slide it into place parallel to the 465 mm extrusions.
- 6. Secure the second assembled side to complete the frame, then verify squareness again with a bubble gauge, adjusting as necessary.
- 7. Then, install the 3D-printed top corner brackets and T-brackets using M5x10 fasteners and M5 T-nuts.
- 8. Add X limit switch mount by attaching limit switch to the mount using M3x16 fasteners and M3 nuts and attaching the mount to the extrusion using M5x10 fasteners and M5 T-nuts.
- 9. Finally, install linear rails by using M3x8 fasteners and M3 T-nuts and following the linear rail alignment procedure outlined in a section below.

6.2 Camera Carriage

- 1. To assemble the Camera Carriage, start with the Enclosure Backplate by attaching the metal plate to the 3D-printed Enclosure Base using four M3x4 fasteners and M3 Hex Nuts.
- 2. Before proceeding further, complete the Belt Assembly system (see below). Next, screw the threaded Camera into the Lens and mount them to the Enclosure Base using a 1/4-20 fastener, placing two 1/4" Washers between the camera and the enclosure.
- 3. Secure the Lens by positioning the Ring Light at the bottom of the Enclosure Base and fastening it with the included M4 Thumb screws, which will later allow for perpendicularity adjustments.
- 4. Mount the Raspberry Pi to the bottom of the Upper Enclosure Cover using four M2.5x16 fasteners and M2.5 Hex Nuts, ensuring its external ports face outward.
- 5. Finally, attach the Upper Enclosure Cover to the Enclosure Base and secure them with the Front Enclosure Cover. Refer to Figure 5 for an exploded view of these steps and Appendix A.2 for a detailed list of fasteners.




Figure 5: Camera carriage assembly; exploded view

6.3 Belt Assembly

The motion in the X-Y plane is controlled using a CoreXY belt configuration. In this motion strategy, the two motors work together to move the carriage in both the X and Y directions. Consequently, the left and right motors are designated as Motor A and Motor B, respectively. Each motor is associated with its own belt, which runs from one side of the carriage assembly, through a set of idlers, around the motor, and back to the other side of the carriage. In this CoreXY application, the A and B loops are stacked vertically in the Z direction, with the B loop always above the A loop. The ends of both loops are clamped between the aluminum backplate of the carriage enclosure and the belt adapter, which has four sets of teeth to mesh with the four belt ends. More information about the CoreXY theory and governing equation can be found at corexy.com.

The assembly of the X-Y belt system is as follows (refer to Appendix A.2 for the detailed fastener list):

- 1. Begin by assembling the Y carriage, which includes the X-axis gantry aluminum extrusion and the left and right brackets.
 - a. Attach the linear rail and trolley to the front face of the aluminum extrusion using M3x8 fasteners and T-nuts.
 - b. Use M5x40 fasteners and hex nuts to attach the upper and lower halves of the brackets around the extrusion.
 - c. Install the idlers, with the front idler consisting of two smooth bearings and the back idler a toothed idler. They are installed using M5x40 fasteners with a washer above and below.
 - d. Attach the Y carriage to the two Y trolleys using four M3x12&x30 fasteners.
- 2. Assemble the XY motor mounts with the stepper motors using M3x8 fasteners and attach the pulleys to the stepper motors.
- 3. Install the motor mounts on the front two vertical extrusions using M5x10 fasteners. Pay attention to the height of the pulleys. Pulley B should be aligned with the upper set of idlers, while Pulley A should be aligned with the lower set of idlers. This can be done roughly at first and then adjusted.
- 4. Assemble idler mounts A and B. Each idler mount contains two toothed idlers and three washers on an M5x40 fastener. The limit switch is installed using M3x16 fasteners.
- 5. Install the idler mounts on the rear left and right vertical extrusions. Idler mount A includes a mounting point for the Y limit switch and should be installed on the left side (closest to Motor A).
- 6. Assemble, but do not install, the belt tensioner. The belt tensioner includes four smooth bearings forming two smooth idlers. An M5x45 fastener runs through the bearings, and a hex nut is used between the bearings to align them.
- 7. Install belt loop A. Belt A runs from the lower left slot on the belt adapter, with the teeth facing the trolley, around the smooth idler on the left Y carriage, around Pulley A, around the lower idler on idler mount A, around the lower idler on idler

- mount B, around the toothed pulley on the right Y carriage, and to the lower right slot on the belt adapter.
- 8. Install belt loop B. Belt loop B follows the opposite path using the upper set of idlers and belt adapter slots.
- 9. Clamp the belts by installing the carriage enclosure. Start by temporarily taping the belt clamp to the X-axis trolley. Then, while holding all four belt ends in position on the belt adapter, use four M3x8 fasteners to attach the carriage enclosure to the trolley and belt adapter, clamping everything together.
- 10. Install the belt tensioner on the back center vertical extrusion using an M5x10 fastener and place the belts in position.

6.4 Platform Assembly




Figure 6: Exploded view of the platform and relevant bill of materials

- 1. First, secure the three Z-axis linear rails to the frame using five M3x8 fasteners and M3 T-Nuts (see Linear Rail Alignment).
- 2. Attach the Z Carriage Brackets to the trolleys on the two front linear rails with four M3x12 fasteners each.
- 3. Next, mount a 370mm aluminum extrusion to the Z Carriage Brackets using two M5x10 fasteners and M5 T-Nuts per side.

- 4. Attach the Z Motor Bracket under the rear linear rail with four M5x10 fasteners and M5 T-Nuts.
- 5. Secure the Z-axis stepper motor beneath the bracket with four M3x12 fasteners, then mount the limit switch using two M3x12 fasteners and M3 Hex Nuts.
- 6. Attach the Ball Screw Bracket to the trolley on the rear linear rail with four M3x12 fasteners. Secure the 255mm aluminum extrusion to the Ball Screw Bracket with six M5x10 fasteners and M5 T-Nuts.
- 7. Fasten the Ball Screw Bracket to the Ball Screw using two M4x14 fasteners with M4 Hex Nuts.
- 8. Join the two aluminum extrusions with corner brackets using four M5x10 fasteners, each with an M5 T-Nut and M5 Washer.
- Attach the Supporting Platform to the extrusions using six M5x10 fasteners, M5
 T-Nuts, and M5 Washers.
- 10. Finally, secure the Breadboard Platform to the Supporting Platform with the M4 fasteners included with the spacers and knobs. (See Appendix A.2 for detailed fastener list)

6.5 Linear Rail Alignment

To install the linear rails, lay the rail onto the extrusion, and use the 3D printed jigs, one on each end, to align the rail properly. The Y and Z rails use the "MGN9H Linear Rail Jig" and the X rail uses the "MGN12H Linear Rail Jig". See the GitHub repository under OpticalModule/CAD/STL Files for the STL files and the source link with pictures.

6.6 Z Motor and Ball Screw Alignment




Figure 7: z-axis assembly top view

Figure: 8 lead screw bracket top view

- 1. First install z-axis assembly onto the aluminum extrusions with 4 M5x10 fasteners and M5 T-nuts.
- 2. Attach one Nema 17 stepper motors with 4 M3x8 fasteners.
- 3. Connect the lead screw and the stepper motor via a coupler.
- 4. Connect the z axis limit switch and secure them with 2 M3x8 fasteners and 2 M3 hex nuts.
- 5. Attach lead screw bracket to linear rail via 4 M3x8 fasteners.
- 6. Attach lead screw flange to lead screw bracket via M3x12 fasteners with M3 hex nuts.
- 7. Adjust the XY position of lead screw with z-axis assembly and lead screw bracket flexible slots until lead screw is perpendicular to stepper motor. Turn the leadscrew until there is little resistance moving from bottom to top.

6.7 Electrical Assembly

6.8 Overview

The below diagram is a representation of the electrical components involved in this project. Since the power supply is rated at 24V and 14.6 A (Max), the voltage needs to be converted to 3 levels to accommodate CNC shield (24V), ring light (12V) and Raspberry pi (5V) respectively. Furthermore, stepper motors and limit switches (NO) are connected to CNC shield via their designated pins. An emergency stop is designed to

terminate and reset the program for CNC shield. Arduino Uno is powered by the Raspberry Pi via a serial connection (USB-A to USB-B). Detail lists of Assembled List and BOM List can be found in Appendix A.6 and A.7.




Figure 9: Electrical components their connections

6.9 E-STOP

The E-Stop is installed upstream of the DC power supply on the live wire between the fused C14 connector and the 120V input terminal of the DC power supply. Ensure ferrules are added to the wires before they are installed in the screw terminals of the E-Stop. After all electrical components are installed, the system should be ground checked before it is powered on to ensure safe operation.

6.10 Buck Converters

The three buck converters power the ring light, the Raspberry Pi, and the CNC shield. The input side of the buck converters should be wired in parallel to the 24V supplied by the DC power supply. For the ring light, the buck converter output should be tuned to 12V using a multimeter and then connected to a DC barrel plug. For the Raspberry Pi, the buck converter should be tuned to 5.1V and then connected to a USB-C breakout

board. For the CNC shield, the buck converter should be tuned to between 15 and 24V and connected directly to the power input terminals of the CNC shield.

7. Software Setup

7.1 Raspberry Pi Settings

Setting Up Raspberry Pi Using Raspberry Pi Imager

Setting up the Raspberry Pi began with using the Raspberry Pi Imager tool. First, download and install Raspberry Pi Imager from the official Raspberry Pi website (https://www.raspberrypi.com/software/). Insert a microSD card or USB drive into your computer and launch the Imager. Click on "Choose OS" and select the preferred Raspberry Pi OS version. For this project, we used the default Raspberry Pi OS (64-bit). Next, click on "Choose Storage" and select the microSD card or USB drive. Next, click on the settings icon to enable SSH, set up a username and password, and configure Wi-Fi if needed. The username and password selected for this project was 'microscope'. Once configured, click "Write" to flash the OS onto the microSD card or USB drive. After completion, eject the storage device and insert it into the Raspberry Pi. Power on the Raspberry Pi, and it will boot into the installed OS.

7.2 Networking

Setting Up Static IP Addresses;

A static IP address was configured on both the Raspberry Pi and Mini PC to establish a direct Ethernet connection between the devices. The Raspberry Pi is physically connected to the mini PC via an Ethernet cable from the Raspberry Pi's Ethernet port to the PC's Ethernet port.

The first step was to configure the computer's static IP address. On Windows, navigate to the Control Panel and open the Network and Sharing Center. Click on Change adapter settings, right-click the Ethernet adapter, and select Properties. Choose Internet Protocol Version 4 (TCP/IPv4), then select 'Use the following IP address'. Enter an IP

address and set the subnet mask to 255.255.255.0. On the Mini PC used for the project, the static IP was configured to be "192.168.1.10". Leave the default gateway blank. Click OK to save the changes.

To configure the Raspberry Pi's static IP address, open a terminal, and type 'ip r | grep default' and save the gateway. Next, type 'nmtui' and select Edit a Connection, followed by Wired Connection. Under IPv4 Configuration, press add and enter the IP address you wish to use. For this project, we used "192.168.1.111" as our address. Also, enter the previously found gateway under Gateway. Afterward, navigate to OK at the bottom of the screen, followed by Back and Quit. Now, type 'sudo systemctl restart NetworkManager' into the terminal. To check if the Static IP address is working, type 'hostname -l' into the terminal.

Once both devices are set up, you can test the connection by opening a terminal or Command Prompt on the computer and running 'ping 192.168.1.111'. A successful connection will return response packets from the Raspberry Pi.

Remote Connection via VNC Viewer:

To remotely access the Raspberry Pi's desktop environment, VNC Viewer was used for this project. To enable the functionality on the Raspberry Pi, first, enable VNC by opening a terminal and running 'sudo raspi-config'. Navigate to Interface Options, select VNC, and enable it. Then, install the RealVNC server if it is not already available using 'sudo apt update' and 'sudo apt install realvnc-vnc-server' in the terminal.

On the PC, download and install VNC Viewer from the official RealVNC website (https://www.realvnc.com/en/connect/download/viewer/). To remotely access the Raspberry Pi, open VNC Viewer and enter the VNC server 'bestcapstone' with the password 'microscope' if the devices are connected to the same wifi network. Otherwise, enter the Raspberry Pi's static IP address (192.168.1.111) in the address bar with the password 'microscope'. If prompted, enter the Raspberry Pi's username, which is also 'microscope'. Once connected, the Raspberry Pi's desktop interface will be

Optical Module Assembly Instructions

Rev. 1, 09-Apr-2025

accessible remotely, allowing full control over the device without needing a monitor, keyboard, or mouse.

7.3 Windows PC

7.3.1 Ethernet

For whichever PC, ensure that it has an ethernet port, or purchase a dongle, so that it can connect to the Raspberry Pi.

7.3.2 Folder Directories

Setup and create folders in an organization like so, and modify gui.py (see <u>GitHub</u>) so that it matches the folder paths you created. This structure is not essential to the program, however, ensure that you are mindful of the folder locations when making modifications to gui.py.

- images
 - o GUI
 - o complete
 - stitching
 - sampling
 - o buffer
 - stitching
 - sampling
 - camera tests

7.3.3 Fiji (ImageJ)


Download Fiji at:

https://imagej.net/software/fiji/downloads


Setting up the Image Stitching Macro in Fiji (ImageJ)

1. Record Macro


Optical Module Assembly Instructions Rev. 1, 09-Apr-2025



2. Image Stitching Plug in



3. Configure Settings (note the Type and Order of scanning do not change)



Directory		Browse
File names for tiles	tile_{ii}.tif	

4. Run the plug-in and save the macro (.ijm file)

The macro written for this project is also available for use at this link:

https://github.com/ManuCapstone1/OpticalModule/tree/main/Programming%20Files/fiji

Appendix

A.1 – Full Parts List

Table 2: List of all parts and tools needed to build the optical module.

Category	Part	Additional Specs or Notes	Website for Specific Parts	Needed Quantity /Size
Electrical	Arduino Uno			1
Electrical	Arduino CNC Shield V3.00	Includes 3 drivers, E-Stop pins		1
Electrical	Mech Endstop V1.2 with Cable			3
Electrical	E-Stop	Twist release, terminal blocks, NC. 24VDC and 10A.		1
Electrical	DC Power Supply	24VDC, ~15A, 50/60Hz		1
Electrical	Power Inlet	250V, 15A		1
Electrical	Raspberry Pi	Pi 5, RAM 2GB		1
Electrical	SD Card for Raspberry Pi	Class 10, 32GB or more		1
Electrical	Stepper Motor	NEMA 17, 1.8deg, 1A		3
Electrical	Stepper Motor Cables	4 Pin, < 60cm long		3
Electrical	Ethernet Cable	90 degree corner facing upwards, cable in direction of locking tab.		1
Electrical	USB C to USB C Cable	>40cm long		1
Electrical	USB C Female Breakout Board PCB	For powering Raspberry Pi from power supply.		1
Electrical	USB-A to USB-B Cable	> 50cm long		1
Electrical	USB A 90 deg corner adapter	For plugging in USB A - USB B cable into raspberry Pi		1
Electrical	microHDMI cable	For displaying Raspberry Pi if needed		1
Electrical	DC-DC Step Down Buck Converter Module	LM2596		3
Electrical	12 AWG Silicone Copper Wires	For wiring power supply to rocker.		>100 cm
Electrical	22 AWG Solid-Core Insulated Wires	For wiring Step Down Converters		> 40 cm
Electrical	Ring Terminal Connectors			10
Electrical	Spade Terminal Connectors			6
Electrical	Cable Mesh Sleeve	>1 m, > 3/8"		1
Electrical	Cable Ties			~5
Electrical	PC	Ethernet Port. Interfaces with Raspberry Pi, storage for images. Runs Windows		1
Electrical	USB Wifi Dongle	Only if needed. To get PC Wifi.		1
Mechanical	Ball Screw and Nut	150 mm long, M8, 2mm per turn	<u>Link</u>	1

Category	Part	Additional Specs or Notes	Website for Specific Parts	Needed Quantity /Size
Mechanical	Lubricant	PTFE Dry lubricant or Superlube		1
Mechanical	Timing Belt	GT2, 6mm		>2 m
Mechanical	Belt Pulley Toothed	GT2, 5mm bore, 6mm, toothed		2
Mechanical	Pulley Idler Toothed	GT2, 5mm bore, 6mm, toothed idler		6
Mechanical	Bearings	F695-2RS, 5x13x4mm, flanged		8
Mechanical	Coupler (motor to balls screw)	5mm to 8mm		1
Mechanical	PLA, PLA +, ABS	Spool (1kg)		2
Mechanical	Linear rails X	MGN12H 300mm		1
Mechanical	Linear rails Y	MGN9H 350mm		2
Mechanical	Linear Rails Z	MGN9H 200mm		3
Mechanical	Spacers and Knobs for Stage	Package of 4 spacers, M4 fastener, and knobs. Typically used for 3D printers		1
Mechanical	Bubble gage			1
Mechanical	Leveling Feet	Adjustable height, max M6 screws, dampening	<u>Link</u>	4
Machined	Breadboard Stage	150x150x13mm	<u>Link</u>	1
Machined	Mounting plate for stage	150x185mm, thickness 4-5mm. See engineering drawing.		1
Machined	Metal plate (camera carriage)	105x80x5mm. See engineering drawing		1
Machined	2020 V-Slot Aluminum Extrusion	370mm x 11 465mm x4 395mm x1 330mm x1 255mm x1		1
Optical	Camera	Raspberry Pi HQ Camera CS 12.3 MP	Link	1
Optical	Lens	100X Industrial Microscope Lens, C/CS-Mount, 0.12X ~ 1.8X	<u>Link</u>	1
Optical	Ring Light	62.5mm inside, 92.5mm outside diameter.	Link	1
Fastener	1/4" Washer			2
Fastener	1/4-20 x 1/2" BHCP or SHCP			1
Fastener	2020 Corner Bracket			34
Fastener	M2.5 Hex Nut			20
Fastener	M2.5x14 BHCP or SHCP			16
Fastener	M2.5x16 FHCP			4
Fastener	M3 Drop T-nut			34
Fastener	M3 Hex Nut			10

Category	Part	Additional Specs or Notes	Website for Specific Parts	Needed Quantity /Size
Fastener	M3x4 BHCP or SHCP			4
Fastener	M3x8 BHCP or SHCP			16
Fastener	M3x8 SHCP			34
Fastener	M3x12 BHCP or SHCP			18
Fastener	M3x16 BHCP or SHCP			4
Fastener	M3x30 BHCP or SHCP			4
Fastener	M4 Hex Nut			2
Fastener	M4x8 BHCP or SHCP			2
Fastener	M4x14 BHCP or SHCP			2
Fastener	M4x20 BHCP or SHCP			2
Fastener	M5 Drop T-Nut			168
Fastener	M5 Hex Nut			12
Fastener	M5 Washers			74
Fastener	M5x10 BHCP			168
Fastener	M5x40 BHCP or SHCP			2
Fastener	M5x40 SHCP			8
Fastener	M5x45 BHCP or SHCP			1
Fastener	M6x14 BHCP or SHCP			2
3D Printed	Enclosure Base V2			1
3D Printed	Enclosure Cover Front			1
3D Printed	Enclosure Cover Upper			1
3D Printed	Belt Adapter			1
3D Printed	Idler Mount B	No limit switch mount.		1
3D Printed	Idler Mount A	Includes limit switch mount.		1
3D Printed	XY Motor Mount	One as is, the other mirrored		2
3D Printed	Z Carriage	One as is, the other mirrored		2
3D Printed	Ball Screw Bracket			1
3D Printed	Z Motor Bracket			1
3D Printed	Housing Base			1
3D Printed	Housing Top			1
3D Printed	Belt Tensioner			1
3D Printed	estop enclosure bottom			1
3D Printed	estop enclosure top			1
3D Printed	Power Supply Housing			1
3D Printed	Sample Holder Pin	Add a spring from a pen to the cylindrical part of the pin.		1
3D Printed	Sample holder clamping plate			1
3D Printed	Sample holder base			1

Optical Module Assembly Instructions Rev. 1, 09-Apr-2025

Category	Part	Additional Specs or Notes	Website for Specific Parts	Needed Quantity /Size
3D Printed	Left Y Carriage Top		GitHub rolohuan	1
3D Printed	Left Y Carriage Bottom		<u>GitHub rolohuan</u>	1
3D Printed	Right Y Carriage Top		GitHub rolohuan	1
3D Printed	Right Y Carriage Bottom		<u>GitHub rolohuan</u>	1
3D Printed	X Limit Switch Mount			1
3D Printed	MGN9 Linear Rail Jig	Print at standard/weaker settings.	<u>Thingiverse</u>	2
3D Printed	MGN12 Linear Rail Jig	Print at standard/weaker settings.	<u>Thingiverse</u>	2
3D Printed	T-bracket		<u>Thingiverse</u>	4
3D Printed	Top Corner Bracket		<u>Thingiverse</u>	4
3D Printed	Ball Screw Knob	Scaled at 97%	<u>Thingiverse</u>	1
3D Printed	Cable Management Clips		<u>Thingiverse</u>	~20
Tools	Metric Allen Keys			1
Tools	Wire Cutter			1
Tools	Crimping Tool			1
Tools	Soldering Iron			1
Machining Tools	Band Saw	If needed, for cutting extrusions.		
Machining Tools	Milling Machine	If needed, for cutting extrusions.		
Machining Tools	Power Drill	If needed, for tapping holes for extrusion leveling feet.		
Machining Tools	Tap and Drill for M6 Holes	If needed, for cutting mounting plate.		
Machining Tools	Waterjet Cutter			

A.2 – List of Fasteners per Part

Table 3 - List of fasteners corresponding to each part.

Row Labels	Sum of # Parts	Sum of Total
Ball Screw Bracket	5	20
M3x8 BHCP or SHCP	1	4
M4 Hex Nut	1	2
M4x14 BHCP or SHCP	1	2
M5 Drop T-Nut	1	6
M5x10 BHCP	1	6
Belt Tensioner	4	5
M5 Drop T-Nut	1	1
M5 Hex Nut	1	2
M5x10 BHCP	1	1
M5x45 BHCP or SHCP	1	1
Enclosure Base	5	15
1/4" Washer	1	2
1/4-20 x 1/2" BHCP or SHCP	1	1
M3 Hex Nut	1	4
M3x4 BHCP or SHCP	1	4
M3x8 BHCP or SHCP	1	4
Enclosure Cover Upper	2	8
M2.5 Hex Nut	1	4
M2.5x16 FHCP	1	4
E-Stop Housing	2	4
M5 Drop T-Nut	1	2
M5x10 BHCP	1	2
Frame	4	238
2020 Corner Bracket	1	34
M5 Drop T-Nut	1	68
M5 Washers	1	68
M5x10 BHCP	1	68
Housing Base	2	32
M2.5 Hex Nut	1	16
M2.5x14 BHCP or SHCP	1	16
Housing Top	2	6
M5 Drop T-Nut	1	3
M5x10 BHCP	1	3
Idler Mount A	6	10
M3 Hex Nut	1	2

Row Labels	Sum of #	Sum of
M3x16 BHCP or SHCP	Parts 1	Total 2
M5 Drop T-Nut	1	2
M5 Hex Nut	1	1
M5x10 BHCP	1	2
M5x40 BHCP or SHCP	1	1
Idler Mount B	4	6
M5 Drop T-Nut	1	2
M5 Hex Nut	1	1
M5x10 BHCP	1	2
M5x40 BHCP or SHCP	1	1
Left Y Carriage	6	18
M3x12 BHCP or SHCP	1	2
M3x30 BHCP or SHCP	1	2
M5 Drop T-Nut	1	3
M5 Hex Nut	1	4
M5x10 BHCP	1	3
M5x40 SHCP	1	4
Linear Rail X	2	10
M3 Drop T-nut	1	5
M3x8 SHCP	1	5
Linear Rail Y	4	28
M3 Drop T-nut	2	14
M3x8 SHCP	2	14
Linear Rail Z	6	30
M3 Drop T-nut	3	15
M3x8 SHCP	3	15
Mounting plate for stage	3	18
M5 Drop T-Nut	1	6
M5 Washers	1	6
M5x10 BHCP	1	6
Power Supply Housing	2	4
M4x8 BHCP or SHCP	1	2
M4x20 BHCP or SHCP	1	2
Right Y Carriage	6	18
M3x12 BHCP or SHCP	1	2
M3x30 BHCP or SHCP	1	2
M5 Drop T-Nut	1	3
M5 Hex Nut	1	4
M5x10 BHCP	1	3

Row Labels	Sum of # Parts	Sum of Total
M5x40 SHCP	1	4
Sample Holder	1	2
M6x14 BHCP or SHCP	1	2
T-Bracket	8	40
M5 Drop T-Nut	4	20
M5x10 BHCP	4	20
Top Corner Bracket	8	72
M5 Drop T-Nut	4	36
M5x10 BHCP	4	36
X Limit Switch Mount	4	8
M3 Hex Nut	1	2
M3x16 BHCP or SHCP	1	2
M5 Drop T-Nut	1	2
M5x10 BHCP	1	2
XY Motor Mount	6	20
M3x8 BHCP or SHCP	2	8
M5 Drop T-Nut	2	6
M5x10 BHCP	2	6
Z Carriage	6	16
M3x12 BHCP or SHCP	2	8
M5 Drop T-Nut	2	4
M5x10 BHCP	2	4
Z Motor Bracket	4	16
M3 Hex Nut	1	2
M3x12 BHCP or SHCP	1	6
M5 Drop T-Nut	1	4
M5x10 BHCP	1	4
Total	NA	644

A.3 – List of Machined Parts

Table 4: List of parts that need to be machined and links to their instructions.

Part	Engineering Drawing	Link to Instructions
Supporting Platform	CAD/Engineering Drawings/supporting platform.pdf	4.2 Plates
Breadboard Platform	CAD/Engineering Drawings/breadboard platform.pdf	4.2 Plates
2020 V-Slot Aluminum Extrusions (Frame)	CAD/Engineering Drawings/430 Main project assy.pdf	4.1 Frame
Enclosure Backplate	CAD/Engineering Drawings/Enclosure Backplate.pdf	4.2 Plates

A.4 – Assembly List - Electronics

Table 5: List of parts that need to be assembled for electronics

Label	Part Type	Properties
24V - 12V	LM2596 step down module	type LM2596 step down
24V - 24V	LM2596 step down module	type LM2596 step down
24V - 5V	LM2596 step down module	type LM2596 step down
A1	Raspberry Pi 5	processor Broadcom 2711, Quad-core Cortex-A76 64-bit SoC @ 2.4GHz; variant Raspberry Pi 5
Arduino CNC Shield V3.00 1	Arduino_CNC_V3.00	variant variant 2
Part1	Arduino Uno (Rev3)	type Arduino UNO (Rev3)
Part2	Adafruit Micro USB Breakout	variant variant 2; part # Adafruit #1833
PI Camera	IC	chip label IC; editable pin labels false; hole size 1.0mm,0.508mm; package DIP (Dual Inline) [THT]; pin spacing 300mil; pins 8; variant variant 6
Power Supply	DC Power Supply	internal resistance 0Ω; max current 14.6A; type dc; voltage 24V
Ring Light	LEDs	package led10mm; variant variant 1
USBTTL1	USB A Plug	variant 1
USBTTL2	USB B Plug	variant 2
X Axis	Stepper Motor - Unipolar	stepper type Unipolar
X limit switch	JL012-13-5.2 Lever (Limit) Switch SPDT + NO + NC	variant mini; part # JL012-13-5.2
Y Axis	Stepper Motor - Unipolar	stepper type Unipolar
Y limit switch	JL012-13-5.2 Lever (Limit) Switch SPDT + NO + NC	variant mini; part # JL012-13-5.2

Optical Module Assembly Instructions Rev. 1, 09-Apr-2025

Z Axis	Stepper Motor - Unipolar	stepper type Unipolar
Z limit switch	JL012-13-5.2 Lever (Limit) Switch SPDT + NO + NC	variant mini; part # JL012-13-5.2

A.5 – BOM - Electronics

Table 8: Bill of Materials for electronics

Part Type	Properties	Amount
LM2596 step down module	Type LM2596 step down	3
Raspberry Pi 5	Processor Broadcom 2711, Quad-core Cortex-A76 64-bit SoC @ 2.4GHz; variant Raspberry Pi 5	1
Arduino_CNC_V3.00	Variant 2	1
Arduino Uno (Rev3)	Type Arduino UNO (Rev3)	1
Adafruit Micro USB Breakout	Variant 2; part # Adafruit #1833	1
Pi Camera	Chip label IC; editable pin labels false; hole size 1.0mm,0.508mm; package DIP (Dual Inline) [THT]; pin spacing 300mil; pins 8; variant 6	1
DC Power Supply	Internal resistance 0Ω; max current 14.6A; type dc; voltage 24V	1
LEDs	Package led10mm; variant variant 1	1
USB A Plug	Variant 1	1
Micro USB B Plug	Variant 2	1
Stepper Motor - Unipolar	Stepper type Unipolar	3
JL012-13-5.2 Lever (Limit) Switch SPDT + NO + NC	Variant mini; part # JL012-13-5.2	3